

## Catalog Number

9512S-0AV-P

**Vials, Screw Top, Glass. Amber, 1.5ml, U Shape. 9-425mm threads, 12x32mm outer dimensions and a fluted stand. For use as an autosampler vial. MICROSOLV Brand. 100/PK.**



### Additional Info

First hydrolytic, Class A, Type N51A, Borosilicate Glass.

Compares to Champagne Style Vials. Features a 15ul residual volume for most autosamplers due to the center draining, inner U shape.

*We recommend to use these vials with [Soft-Guard™](#) 9-425mm Screw Caps.*

**Price**  
**\$91.83**

## Additional Info

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$91.83**