

Catalog Number

9512S-0CV-T

Vials, Screw Top, Glass. Clear, 1.2ml, MRQ™, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. AQ Brand. 100/PK.



Additional Info

MRQ Vials: Optimal Sample Visibility and Recovery

Our MRQ vials are designed for optimal sample visibility and recovery. Made from first hydrolytic, Class A, Type-33 Borosilicate Glass, these vials have a fill volume of 1.2ml and an exceptionally low residual volume of 2 μ l. This is achieved through the center-draining effects of smooth walls (no ledges) and a conical bottom with a well slightly larger than autosampler vial needles. This design ensures that bottom entry port needles can recover the maximum amount of sample during use.

Price

\$204.95

Additional Info

The bottom depths of the wells are meticulously designed and quality-controlled to prevent needles from hitting the bottom of the vials, ensuring compatibility with the most critical autosamplers.

Certificate of Conformance is available upon request.

Click [HERE](#) for a schematic drawing of the vial with specifications. This product is manufactured in an ISO 9001:2015 Certified facility.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$204.95