

## Catalog Number

9512S-3DCP-PX

**Vial & Cap Kit.** 100 EA 1.8ml, Screw Top, Hydrophobic, Clear Autosampler Vials and matching Black Caps with Clear Silicone Rubber / PTFE Pre-Slit Septa. RSA-Pro X Brand. 1 PK.



### Additional Info

These 1.8ml RSA, Glass, Max Recovery Autosampler Vial & Cap Kit Packs (*Easy Purchase Packs*) are supplied with [Clear Center Draining RSA Vials](#) that have been treated with [RSA-Pro X™ Technology](#) to prevent their already pristine surface from adsorbing stubbornly sticky compounds that have affinity for the glass structure. The Hydrophobic surface is not only non-adsorbing for most compounds but is extremely stable even in the presence of Water. The hydrolytic stability will last for years compared to Silanized Vials that will usually degrade in the presence of Water. This product is manufactured in an ISO 9001:2015 Certified facility.

The AQR™ Brand, matching [Black Screw Caps](#) have the most inert Clear, Silicone Rubber / Clear, PTFE fitted Pre-Slit Septa available; for total Vial and Cap inertness. This Kit is often chosen for proteins, peptides, enzymes, natural products & unknowns especially when working in parts per billion with LCMS or HPLC.

*These treated vials are only sold in RSA-Pro X™ Kits. Click [HERE](#) for replacement Caps.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$379.60**

Printed On 2024/05/07