

Catalog Number

95296-1mL-3P

Micro-Sample Management System. This set includes 96, 1000ul Precision, Glass Inserts with wide conical points in a Rack with a Base and a Pre-Slit, Silicone Rubber, Plug Style Sealing Mat. U-2D Brand.



Additional Info

This Set is a convenient way to purchase all that is needed to manage your Micro-Samples or use this set in a 96 Well Plate configuration. 1ml Inserts are 1st Hydrolytic, Class A Type 33 Borosilicate Glass. Use with Aqueous Solutions.

These 1ml Inserts have a wide conical point for fast sample filling to minimize bubble formation and has outer dimensions of 7.72mm x 44.85mm.

Set includes [Pre-Racked Inserts](#), [Base](#) & [Pre-Slit Sealing Mat](#).

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$137.23

Printed On 2024/05/06