


Catalog Number

95296-1mL-P-CS

Inserts, RSA-Pro X, in Racks. Includes 1,152 clear, 1mL surface treated glass Inserts with wide conical Points in 12 Racks of 96. U-2D Brand. 1 CS.



Additional Info

These RSA-Pro X™ treated Inserts are supplied in a Rack for easy filling and managing. [The Inserts are 1mL precision bottom Glass Inserts](#) that are treated with [RSA-Pro X™ Surface Technology](#) to upgrade the pristine surface to further resist adsorbing stubbornly sticky compounds that have affinity for the glass structure. The Hydrophobic glass surface is not only non-adsorbing for most compounds but is extremely stable even in the presence of Water. The hydrolytic stability will last for years compared to silanized vials that will degrade in the presence of Water. This product is manufactured in an ISO 9001:2015 Certified facility.

Use single or multi channel pipettors to fill the Inserts in the [U-2D™ Rack](#) which can be placed in a matching Base. Inserts can be covered and sealed with a Solvent Resistant Mat (organic

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,605.04

Printed On 2024/08/06

Additional Info

solvents) or a Silicone Rubber (aqueous) Mat and used for storage or directly in an autosampler that accepts, standard ANSI format, 96 Well Plates.

RSA-Pro X™ is often needed to prevent standards, proteins, peptides, enzymes, natural products and other precious compounds from sticking to borosilicate glass. U-2D™ Base and Mats are sold separately.

These 1ml Inserts have a wide, conical point for fast filling and to minimize bubble formation. Inserts outer dimensions: 7.72mm x 44.85mm.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$1,605.04

Printed On 2024/08/06