

## Catalog Number

95296-1mL-PK

**Inserts, in Racks.** Includes 96, Clear, 1000ul Glass Low Volume Precision Inserts with wide conical points. For use in the Micro-Sample Management System. U-2D Brand. 6/PK.



### Additional Info

Inserts are 1st Hydrolytic, Class A, TYPE 33 Borosilicate Glass made for use with the U-2D Micro-Sample Management System.

These 1ml Inserts have a wide conical point for fast sample filling to minimize bubble formation and have outer dimensions of 7.72mm x 44.85mm.

Preloaded in Racks. Click [HERE](#) for case quantity. *Bases and Sealing Film or Mats sold separately.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$760.71**

Printed On 2024/05/06