

## Catalog Number

95296-30R-CS

**Mats, Sealing, Chemical Resistant Polymer, Pre-Slit, Plug Style for 96, 350ul Inserts. U-2D Brand. 5/CS.**



### Additional Info

Made of a chemically resistant polymer, this Mat can be used with typical HPLC organic solvents. Designed and molded to seal up to 96 each 350ul Inserts in the U-2D Micro-Sample Management System. This mat is Pre-Slit and has alphanumeric markings on the top side of the mats, along the edges for identification of the Inserts.

Click [HERE](#) for individual quantity.

**Price**  
**\$132.13**

## Additional Info

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$132.13**