

## Catalog Number

95296-350-1R

**Micro-Sample Management System.** This set includes 96, 350ul Precision, Glass Inserts in a Rack with a Base and a Solvent Resistant, Polymeric Plug Style Sealing Mat. U-2D Brand. 1/PK.



### Additional Info

This Set is a convenient way to purchase all that is needed to manage your Micro-Samples or use this set in a 96 Well Plate configuration. Inserts are 1st Hydrolytic, Class A Type 33 Borosilicate Glass. Use this set with most HPLC compatible, organic solvents.

These 350ul inserts can easily be removed from the rack and placed in a 12x32mm, 2ml AQ™ Brand autosampler vial for use in GC or HPLC.

Set includes [Pre-Racked Inserts](#), [Base](#) & [Non-Slit Sealing Mat](#).

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$127.94**

Printed On 2024/05/02