

Catalog Number

95296-350-3R

Micro-Sample Management System. This set includes 96, 350ul Precision, Glass Inserts in a Rack with a Base and a Pre-Slit Solvent Resistant, Polymeric Plug Style Sealing Mat. U-2D Brand. 1/PK.



Additional Info

This Set is a convenient way to purchase all that is needed to manage your Micro-Samples or use this set in a 96 Well Plate configuration. Inserts are 1st Hydrolytic, Class A Type 33 Borosilicate Glass. Use this set with most HPLC compatible, organic solvents.

These 350ul inserts can easily be removed from the rack and placed in a 12x32mm, 2ml AQ™ Brand autosampler vial for use in GC or HPLC.

Set includes [Pre-Racked Inserts](#), [Base](#) & [Pre-Slit Sealing Mat](#).

Price

\$108.03

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$108.03