

Catalog Number

95296-500-3R

Micro-Sample Management System. This set includes 96, 500ul Precision, Glass Inserts in a Rack with a Base and a Pre-Slit Solvent Resistant, Polymeric Plug Style Sealing Mat. U-2D Brand. 1/PK.



Additional Info

This Set is a convenient way to purchase all that is needed to manage your Micro-Samples or use this set in a 96 Well Plate configuration. Inserts are 1st Hydrolytic, Class A Type 33 Borosilicate Glass. Use this set with most HPLC compatible, organic solvents.

Set includes [Pre-Racked Inserts](#), [Base](#) & [Pre-Slit Sealing Mat](#).

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$143.85

Printed On 2024/04/29