

Catalog Number

95296-I-500

Inserts, Low Volume, Glass. Clear, 500ul with a wide, conical point for use in the Micro-Sample Management System. Replacement Inserts. U-2D Brand. 100/PK.



Additional Info

Inserts are 1st Hydrolytic, Class A, TYPE 33 Borosilicate Glass made for use with the U-2D Micro-Sample Management System Rack and have outer dimensions of 7.72mm x 31.30mm

Also available pre-loaded in U-2D Brand Racks of 96. Click [HERE](#) for more information.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$114.23

Printed On 2024/05/03