

Catalog Number

95296-LG-10P-CS

Mats, Sealing, Large, Silicone Rubber, Plug Style for 96, 500 and 1000ul Inserts. U-2D Brand. 5/CS.



Additional Info

Mats are made of ultra pure Silicone Rubber, designed and molded to seal up to 96 each 500ul and 1ml Inserts in the U-2D Micro-Sample Management System. These mats are not Pre-Slit and have alphanumeric markings on the top side of the mats, along the edges for identification of the Inserts. Use with Aqueous Solutions.

Click [HERE](#) for single unit quantity.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$103.74

Printed On 2024/06/20