

## Catalog Number

95296-LG-30P-CS

**Mats, Sealing, Large, Silicone Rubber, Pre-Slit, Plug Style for 96, 500 and 1000ul Inserts. U-2D Brand. 5/CS.**



### Additional Info

Mats are made of ultra pure Silicone Rubber, designed and molded to seal up to 96 each 500ul and 1ml Inserts in the U-2D Micro-Sample Management System. These mats are Pre-Slit and have alphanumeric markings on the top side of the mats, along the edges for identification of the Inserts. Use with Aqueous Solutions.

Click [HERE](#) for single unit quantity.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$129.91**

Printed On 2024/05/08