

Catalog Number

95296-LG-30R-CS

Mats, Sealing, Large, Chemical Resistant Polymer, Pre-Slit, Plug Style
for 96, 500 or 1000ul Inserts. U-2D Brand. 5/CS.



Additional Info

Mats are made of a chemically resistant polymer that can be used with typical HPLC organic solvents. Designed and molded to seal up to 96 each 500ul or 1ml Inserts in the U-2D Micro-Sample Management System. These mats are Pre-Slit and have alphanumeric markings on the top side of the mats, along the edges for identification of the Inserts.

Click [HERE](#) for single unit quantity.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$129.23

Printed On 2024/05/03