

Catalog Number

95296-MARK-FS

Marker, permanent. Extra-Fine tip, 0.25mm. Silver Ink for plasticware. 1 EA.



Additional Info

Mark in bold, opaque silver color on projects and more with this marker. The Extra Fine tip features a 0.7 mm pin-type plastic tip and

Valve-action with rich, smooth-flowing paint.

Archival quality, waterproof, and fade resistant on most surfaces including, plastic, wood, porcelain, paper and metal. Permanent adhesion on most surfaces- except for glass and ceramics.

Low-odor and xylene free. Meets ASTM and ACMI non-toxicity standards.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$11.08

Printed On 2024/07/27