

## Catalog Number

95296-MARK-FS

**Marker**, permanent. Extra-Fine tip, 0.25mm. Silver Ink for plasticware. 1 EA.



### Additional Info

Mark in bold, opaque silver color on projects and more with this marker. The Extra Fine tip features a 0.7 mm pin-type plastic tip and

Valve-action with rich, smooth-flowing paint.

Archival quality, waterproof, and fade resistant on most surfaces including, plastic, wood, porcelain, paper and metal. Permanent adhesion on most surfaces- except for glass and ceramics.

Low-odor and xylene free. Meets ASTM and ACMI non-toxicity standards.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$11.08**

Printed On 2024/05/24