

Catalog Number

95296-MARK-OB-P

Marker, Laboratory, Semi - Permanent. 1.5mm Line. 13ml, Blue Ink. 4/PK.



Additional Info

Mark in bold, blue color on projects and more with this marker. This indelible, oil-based ink works on all surfaces. Once the ink is dried it presents a semi-permanent finished on Hydrophobic Glass Vials and Inserts. Click [HERE](#) for Hydrophobic RSA-Pro X™ Glass Vials and Inserts.

Laboratory Markers Features:

- Tip Type: Fine Point
- Marking Width: 1.5mm
- Marking Temperature Range: -30° to 60° C
- Dry Time: 10 Seconds

Directions for use: Shake vigorously with cap on. Press the tip on a surface that can be marked until inked. Replace the cap after use.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$47.31

Printed On 2024/05/06