

## Catalog Number

95296-U2D

**Inserts, in Racks with a Base.** Includes 96, Clear, 350ul Glass Low Volume Precision Inserts with wide conical points. For use in the Micro-Sample Management System. U-2D Brand. 1 EA.



### Additional Info

Inserts are 1st Hydrolytic, Class A, TYPE 33 Borosilicate Glass made for use with the U-2D Micro-Sample Management System or in 2ml, 12x32 AQ Brand Autosampler Vials, Pre-Loaded in a Rack with a Base.

These Inserts have a wide conical point for fast sample filling and minimizes bubble formation. *Sealing Film or Mats sold separately.*

**Price**

**\$118.15**

## Additional Info

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$118.15**