

## Catalog Number

9532C-5DCP

**Vial & Cap Kit** includes 100 EA 300ul, Fused Insert, Snap Top, Clear Deactivated Autosampler Vials with Write on Patch and fill lines and 100 EA matching Snap Caps with ultra pure, PTFE, fitted in the Blue Caps.



### Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps. The kit packaging will open and re-close making it useful for keeping vial and cap types matched, organized and clean during use.

Click on links below for more Specifications for the 11mm Vials and Caps:

Supplied with Silanized / Deactivated Vials: Cat. No [9532C-0DCV](#) and Caps: Cat No. [9509C-50CB](#)

Price

\$291.12

## Additional Info

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$291.12**