

## Catalog Number

## 9532S-0PV-K

**Vials, Screw Top, Plastic.** Clear, 300ul MicroVial, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial. MicroSolv Brand. 1000/CS.



### Additional Info

This narrow mouth, vial is manufactured with virgin, medical grade polypropylene and has a clear, translucent, built-in insert with a 300ul total fill capacity. Residual volume can be as low as 2ul depending on the autosampler and viscosity of samples used.

Popular choice for samples that react with glass and often used in bio-analysis with aqueous samples.

Click [HERE](#) for smaller pack quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$236.66**

Printed On 2024/06/20