

## Catalog Number

## 9532S-1ECP-WTRS

Vial & Cap Kit Includes 100 EA 300ul, Fused Insert with a wide tip, Screw Top, Clear RSA Autosampler Vials with Write on Patch and fill lines and 100 EA matching Screw Caps with White Silicone Rubber / Red PTFE



### Additional Info

Easy Purchase Pack™ Vial and Cap Kits are a convenient way to order and keep inventory of autosampler vials with their matching caps organized. The kit packaging will open and re-close making it useful for keeping vial and cap types together and clean during use and storage.

Click on links for more Specifications for the 9mm Vials and Caps: Supplied with Reduced Surface Activity Vials (RSA™), Cat. No. [9532S-WCV-WT-RS](#) and Soft-Guard™ Brand Caps, Cat No. [9502S-E-10C-B](#) .

**Price**

**\$225.58**

## Additional Info

*This is one of our Most Popular Vial Kits.*

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$225.58**