

## Catalog Number

9532S-E-0CV

**Vials, Screw Top, Glass. Clear / Clear Insert, 300ul Fused Insert™, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial. BASIK Brand. 100/PK.**



### Additional Info

Vials and inserts are first hydrolytic, Class A, Borosilicate Glass, with a wide bottom, 300ul, inserts fused to the bottom of the vial. The wide bottom helps prevent bubble formation when filling quickly or with automation.

Our most economical Fused-Insert Vials. For budget minded laboratories.

**Price**  
**\$146.95**

## Additional Info

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$146.95**