

## Catalog Number

9532S-MS

**Vials, Screw Top, Plastic.** Clear, 300ul MicroVial, LCMS Compatible, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial. AQ Brand. 100/PK.



### Additional Info

This narrow mouth vial is manufactured with ultra-pure, virgin, medical grade polypropylene and is translucent with a built-in insert that has 300ul total fill capacity. Residual volume can be as low as 2ul depending on the autosampler and viscosity of samples used.

This specially formulated plastic is used in our special molding process where nothing is added to the process producing a very unique, plastic vial is that compatible with LCMS (extractables free in most LCMS solvents). This product is manufactured in an ISO 9001:2015 Certified facility.

**Price**  
**\$46.89**

## Additional Info

Popular for samples that react with glass and often used in bio-analysis with aqueous samples.

Click [HERE](#) for case quantities.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$46.89**