

## Catalog Number

9534S-3DCP

**Vial & Cap Kit** Includes 100 EA 3.8ml, Maximum Recovery, Screw Top, Clear Deactivated Library Storage Vials and 100 EA matching Screw Caps with Silicone Rubber / PTFE, Pre-Slit Septa bonded in the Black



### Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps.

Click on links below for more Specifications for the 13mm Vials and Caps:

Supplied with Silanized / Deactivated Vials: Cat. No [9534S-0DCV](#) and Caps: Cat No. [9504S-30C-BK](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$383.69**

Printed On 2024/05/05