

Catalog Number

9534S-E-0CV-K

Vials, Screw Top, Glass. Clear, 3.8ml, Max Recovery™, 13-425mm threads, 14.75x45mm outer dimensions for 13mm screw caps. For use as an autosampler vial or storage container. BASIK Brand. 1000/CS.



Additional Info

First hydrolytic, Class A, pharmaceutical grade, Borosilicate Glass. This 3.8ml vial has the low residual volume of 18ul.

Low priced for budget minded laboratories.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$2,703.99

Printed On 2024/04/27