

Catalog Number

9534S-SPA13

Vial & Cap Kit includes 100 EA 3.8mL, Max Recovery™, Screw Top, Amber Autosampler Vials and 100 EA matching Screw Caps with a Silicone Rubber / PTFE Liner in the Black solid Caps. AQ Brand Easy



Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps. The kit packaging will open and re-close making it useful for keeping vial and cap types matched, organized and clean during use.

Click on links below for more Specifications for the 9mm Vials and Caps:

Supplied with Vials: Cat. No [9534S-0AV](#) and Caps: Cat No [98100-SP13-K](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$399.88

Printed On 2024/05/04