

## Catalog Number

9572S-PP

**Vials, Screw Top, Plastic.** 700ul, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial. MICROSOLV Brand. 100/PK.



### Additional Info

These injection molded, translucent, plastic vials are made with virgin polypropylene and chosen for bio-analytical reasons.

Click [HERE](#) for case quantities.

**Price**  
**\$30.94**

## Additional Info

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$30.94**