

## Catalog Number

9582S-02N

**Inserts, Low Volume** for Vials, Glass. Clear, 200ul fill volume with a conical point and 5x29mm outer dimensions with plastic springs attached. For use in 8mm, 2ml autosampler vials. MicroSolv Brand.



### Additional Info

First hydrolytic, Class A, Type N51A, Borosilicate Glass. Features a precision point conical tip for excellent recovery of sample. Residual volume is <1ul with most autosamplers with a 100ul useable volume. Polypropylene springs are attached to keep the inserts center aligned and reduce chances of the needle breaking the inserts during use.

Use with 8mm, narrow mouth Screw Top or narrow mouth Serum Finish Crimp Top 12x32mm autosampler vials.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$128.83**

Printed On 2024/05/24