

Catalog Number

9582S-05N-K

Inserts, Low Volume for Vials, Glass. Clear, 300ul fill volume with a flat bottom and 5x31mm outer dimensions. For use in 2ml autosampler vials. MicroSolv Brand. (10 packs of 100) 1,000/CS.



Additional Info

First hydrolytic, Class A, Type N51A, Borosilicate Glass. Features a flat bottom. Residual volume is 8ul with most autosamplers.

Use with standard opening, 8mm Screw Top or 11mm Serum Finish Crimp Top 2ml, 12x32mm autosampler vials. *Lower cost for budget minded laboratories.*

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$322.58

Printed On 2024/05/18