

## Catalog Number

9582S-0AV

**Vials, Screw Top, Glass. Amber, 2ml, 8-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. BASIK Brand. 100/PK.**



### Additional Info

First hydrolytic, Class A, Borosilicate Glass, amberized. Features a narrow mouth which may be useful to reduce evaporation and requires an 8mm screw cap.

Amber vials are often used for light sensitive compounds and considered low actinic.

*Lowest cost in class.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$33.76**

Printed On 2024/04/29