

## Catalog Number

9582S-20C-K-M

**Caps, Screw Top**, with PTFE / Silicone Rubber / PTFE Septa, fitted in "knurled", polypropylene, Black Caps.. For use with 8-425mm thread, 12x32mm autosampler vials. BASIK Brand. 1000/CS.



### Additional Info

This septa was designed for specific, older instruments such as Agilent 1050 and to provide extra cleaning when "wiping" a needle. Also provides a small degree of better resealing after injection compared to Silicone Rubber / PTFE Septa.

Click [HERE](#) for smaller pack quantities.

**Price**  
**\$298.83**

## Additional Info

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$298.83**