


Catalog Number

96008-10-5010

Septa, 8mm, White Silicone Rubber / Beige PTFE that is 0.060 inch thick and 0.321 inch diameter for use in 8-425 mm thread caps. Caps sold separately. MICROSOLV Brand. 1000/CS.



Additional Info

Used for both our HPLC or GC, 8mm screw caps.

Price
\$228.98

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$228.98