


## Catalog Number

96010-10-1501

**Septa, 10mm**, White Silicone Rubber / Beige PTFE that is 1.5mm thick for use in 10-425mm thread caps. Caps sold separately. MicroSolv Brand. 1000/CS.



### Additional Info

Use for HPLC or GC.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$645.84**

Printed On 2024/05/07