




Catalog Number

96020-1X-9010

Septa, 20mm, Headspace, Clear Silicone Rubber / White PTFE. For use in 20mm Crimp Caps, sold separately. BASIK Brand. 1000/CS.



Additional Info

Silicone Rubber / PTFE Septa 100mil, 90mil rubber and 10mil PTFE. This Septa is good for higher temperature analyses.

Temperature range from -60° to 200°C.

Use in Aluminum or Steel Crimp Caps for Headspace autosamplers. Supplied in Polyethylene Bags.

Price

\$763.71



Additional Info

MICROSOLV

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$763.71