

## Catalog Number

97004-01

**Vial & Cap Kit** Includes 144 EA 4ml, Screw Top, Clear Storage Vials and 144 EA matching, Attached Screw Caps with a Rubber / PTFE Liner fitted in the Black, Knurled, Phenolic Caps. BASIK Brand. Easy Purchase



### Additional Info

These 4ml Vials are made with 1st Hydrolytic, TYPE 33 Borosilicate Glass with 13-425mm threads that a supplied in egg crate style corrugated tray, boxes.

Closed, or solid, Phenolic caps are supplied and can tolerate higher temperatures than polypropylene and are harder.

*Caps and Vials are supplied assembled.*

**Price**

**\$493.16**

## Additional Info

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$493.16**