

Catalog Number

97010-CV05

Vials, Push Top, Storage, Glass, Clear, 5ml. Economical, with ND18 tops. For use with Econo-Storage Push Top Caps. BASIK Brand. 100/PK.



Additional Info

These Vials are made with 3rd Hydrolytic Glass and quickly open and reseal due to the "Snap" design. These 5ml Vials have a total fill volume of 9.1ml with 40x20mm outer dimensions. Not made for Crimp Caps.

Use [97010-SC18](#) Caps with these Vials, sold separately. These vials are not intended for Analytical Testing and should only be used for Storage.

Low cost for budget minded laboratories.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$32.91

Printed On 2024/04/29