

Catalog Number

97015-AV-08

Vials, Screw Top, storage, glass, amber, 8ml, with 15-425mm screw threads. For use as a storage vial. MicroSolv brand. 100/PK.



Additional Info

These Storage Vials have 61x16.6mm outer dimensions and a total fill volume of 8.9ml and can be used for long term storage with a proper cap for light sensitive compounds.

Manufactured with amberized, low actinic First Hydrolytic, Class A, Type N51A, pharmaceutical grade borosilicate glass. Use with caps with 15-425 or ND15 screw threads which are sold separately.

Click [HERE](#) for clear vials.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$71.18

Printed On 2024/04/29