

Catalog Number

97015-AV-12

Vials, Screw Top, Storage, Glass, Amber, 12ml, with 15-425mm Screw Threads. For use as a Storage Vial. MicroSolv Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type N51A, Pharmaceutical Grade, Borosilicate Glass, amberized. Amber vials are often used for light sensitive compounds and are considered low actinic.

These Storage Vials have 66x18.5mm outer dimensions and a total fill volume of 12ml and can be used for long term storage with a proper cap.

Use with Caps that have 15-425 or ND15 screw threads, sold separately.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$72.50

Printed On 2024/08/06

Additional Info

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$72.50

Printed On 2024/08/06