

## Catalog Number

## 97015-CV-08

**Vials, Screw Top, storage, glass, clear, 8ml, with 15-425mm screw threads. For use as a storage vial. MicroSolv brand. 100/PK.**



### Additional Info

These Storage Vials have 61x16.6mm outer dimensions and a total fill volume of 8.9ml and can be used for long term storage with a proper cap.

Manufactured with First Hydrolytic, Class A, Type N51A, pharmaceutical grade, borosilicate glass. Use with caps with 15-425 or ND15 screw threads which are sold separately.

Click [HERE](#) for amber vials.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$67.89**

Printed On 2024/02/25