

## Catalog Number

97015-CV-12

**Vials, Screw Top, Storage, Glass, Clear, 12ml, with 15-425mm Screw Threads. For use as a Storage Vial. MicroSolv Brand. 100/PK.**



### Additional Info

First hydrolytic, Class A, Type N51A, Pharmaceutical Grade, Borosilicate Glass.

These Storage Vials have 66x18.5mm outer dimensions and a total fill volume of 12ml and can be used for long term storage with a proper cap.

Use with Caps that have 15-425 or ND15 screw threads, sold separately.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$66.39**

Printed On 2024/05/24