

## Catalog Number

97018-CV-16

**Vials, Screw Top, Storage, Glass, Clear, 16ml, with 18-400mm Screw Threads.** For use with For use as a Storage Vial. MicroSolv Brand.  
100/PK.



### Additional Info

First hydrolytic, Class A, Type N51A, Pharmaceutical Grade, Borosilicate Glass.

These Storage Vials have 71x20.6mm outer dimensions and a total fill volume of 17.4ml and can be used for long term storage with a proper cap.

Use with Caps that have 18-400 or ND18 screw threads, sold separately.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$89.18**

Printed On 2024/05/07