


Catalog Number

97041-11

Septa, 22mm, Natural Silicone Rubber / White PTFE. Ultra Low Bleed.
For use in 24mm Caps, sold separately. BASIK Brand. 1000/CS.



Additional Info

Silicone Rubber / PTFE Septa 2.5mm thick. This Septa is good for storage or EPA Testing due to the high purity of the rubber and PTFE. Temperature range from -60° to 200°C.

Use in screw caps that have 24mm or ND2 threads.. Supplied in a polyethylene Bags.

Price

\$1,021.00

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,021.00