


Catalog Number

97041-SP24

Caps, Screw Top, Solid, with Natural Silicone Rubber / Natural PTFE Liner, fitted in knurled polypropylene, White Caps. For use with 24-414mm threaded EPA Style vials. ND24 threads. BASIK Brand. 100/PK.



Additional Info

Liner is 3.18mm (0.125") thick and can be used for long term storage with any 30, 40 or 60 EPA Style vial that has 40-414 threads.

Solid top cap does not have an opening and is designed for storage or transportation.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$115.94

Printed On 2024/07/27