


Catalog Number

98100-SP20-W

Caps, Screw Top, Solid, with Silicone Rubber / PTFE Liner, fitted in knurled polypropylene, White Caps. For use with ND20 thread, storage Vials or Bottles. BASIK Brand. 100/PK.



Additional Info

Silicone Rubber / PTFE Liner can be used for use with Storage Vials / Bottles with 20-400mm screw threads. [97020-CV-20](#).

Also, used with heavy wall reaction vials.

Price
\$77.13

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$77.13