

## Catalog Number

98100-SP20-W

**Caps, Screw Top, Solid, with Silicone Rubber / PTFE Liner, fitted in knurled polypropylene, White Caps. For use with ND20 thread, storage Vials or Bottles. BASIK Brand. 100/PK.**



### Additional Info

Silicone Rubber / PTFE Liner can be used for use with Storage Vials / Bottles with 20-400mm screw threads. [97020-CV-20](#).

Also, used with heavy wall reaction vials.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$77.13**

Printed On 2024/05/02