

Catalog Number

MI10-L9010

LABeler Lab Printer, portable, Tarsus Brand. 1 EA.



Additional Info

The LABeler™ label maker eliminates the need for a lab pen. This portable device can be used to type in the desired copy and print within seconds. The label can be cut by pressing the cut lever and is ready for application. The LABeler is ideal for label identification of any laboratory vessel for nearly any laboratory application.

Standard label tapes (6, 12 or 24mm wide labels) can be autoclaved, boiled or even frozen to temperatures as low as -50°C. [Special Ultra-Low](#) (U) label tapes are also available for storage in -80°C environments or for liquid nitrogen applications.

Supplied with 115V AC/DC power supply with 9V output but can run on 6 AA batteries as well.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$243.75

Printed On 2024/08/06

Additional Info

US plug supplied.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$243.75

Printed On 2024/08/06