

Catalog Number

PQA-GVS-01

Gradient Visualization Test Mix Solution. Replacement. 30 ml.
Chemical Solutions brand. 1 EA.



Additional Info

This NIST traceable replacement solution is for testing the precision and accuracy of your HPLC gradients. Use with the Chemical Solutions PQ or HSQ kits. Supplied in a 30 ml plastic bottle and cap.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

Price
\$291.44

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$291.44